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## IMPORTANT INFORMATION FOR THE USE OF THE KONE-LOG IMPLANT

### DESCRIPTION

The subtalar implant Kone-Log is a one-piece implant for treating hyperpronation foot made of titanium alloy (Ti-6AL-4V ELI). Kone-Log is available in 6 different sizes with a diameter from 6 to 13 mm in diameter. The implants are available in sterile single package. The instruments to facilitate the placement of the implants are packed in a autoclavable box. The box includes : series of cannulated test calibrators, the guide wire, the cannulated screwdriver for insertion and the cannulated screw driver for extraction. The instruments are made of stainless steel and polymer material.

### CONDITIONS OF USE

The implant must be used by practitioners surgeons who have received the necessary information about the implant of this device.

For correct use of the Kone-Log implant the surgeon should be aware that **the correct selection and sizing of the implant is crucial**. The choice of the size, shape and structure of the implant increases the chances of success. The implants should be placed with care and adequate bone stock.

### INDICATIONS

The subtalar Kone-Log implant is indicated for use for the purpose of surgical treatment of hyperpronation foot and subtalar joint stabilization in pediatric patients, because no two foot measurements are comparable, research should consider developing recommendations on pediatric foot measurement using valid and reliable assessment tools. . This system is designed to block the movement of the heel forward, downward or in a medial direction to obtain the restoration of normal subtalar joint mobility while preventing the excessive pronation and the consequences.

In the indications for use include the following:

- pronounced pronation;
- calcaneal valgus deformity;
- plantar flexion of the heel;
- not corrected with orthotic treatment in the long term;
- congenital flat foot and painful;
- repair of coalescence tarsal;
- subtalar instability;
- failure of the posterior tibial tendon;
- paralytic flat foot.

### CONTRAINDICATIONS

The use of subtalar implant Kone-Log is contraindicated in the following cases:

- Active or suspected infection or in immunocompromised patients;
- In patients with previously sensitized to titanium;
- In patients with inadequate bone stock;
- In patients with superstructural alignment deformities or in patients with certain metabolic disorders;
- patient's disagreement.

The evaluation of the contraindications is however postponed to the competence of the surgeon, who decides whether they are to be considered relative or absolute.

### PERFORMANCES

The mechanical characteristics of Kone-log implant must allow hyperpronation correction of foot and subtalar joint stabilization. Performance may be hindered by a bad use of the material or an inadequate post-operative patient monitoring.

### SIDE EFFECTS

The patient should be put to the knowledge by the surgeon general risks associated with surgery before the operation itself.

- Current infection;
- Pain, discomfort or abnormal sensations due to presence of the implant;
- Metal sensitivity or allergic reaction to a foreign body;
- Dislocation of the implant;
- Loosening of the implant;
- Alignment belated correction;
- Reduction in bone density due to stress shielding;
- Bursitis.

### PRECAUTIONS

The Kone-Log implants are supplied STERILE. Sterilization by Ethylene Oxide (EO) in compliance with UNI EN ISO 11135-1. Use only if the package is intact and undamaged. Do not re-sterilize.

### SINGLE USE.

### Preoperative

A good knowledge of the implant and the surgical technique are essential. The usage range of the device is entrusted to the doctor competent surgeon according to the indications and contraindications and the patient's anatomy.

### Intra-operative

The instructions for use of Kone-Log implant must be followed carefully on the surgical technique. It 'not recommended to associate our implants to other similar devices.

### Postoperative

A number of precautions must be given to the patient by the surgeon: limitation of his physical activities, compliance with the instructions of the surgeon.

### CONTROL

Control of assembly and through visual image intensifier; Clinical correction control.

### Nuclear Magnetic Resonance (NMR)

Security and compatibility within the MR environment with the devices described in this instruction for use have not been evaluated. The implant described in this instruction for use have not been tested with regard to the heating or migration in the MR environment.

### PACKAGING

Each Kone-Log implant is provided STERILE and SINGLE USE: will be used directly in the operating room. Reuse may cause a re-infection or cross-infection, or impair the functional performance of the device. Sterility is guaranteed only if the packaging is intact and undamaged. Do not re-sterilize.

### ASSEMBLY OF THE DEVICE

The use of the Kone-Log implant is made by the surgeon following the surgical technique, recommended on the surgical technique.

### RADIOACTIVE DEVICES

The Kone-Log implant is not radioactive.

### CHANGE PERFORMANCES

This implant is expected to perform the treating hyperpronation of the foot and of the subtalar joint stabilization. A healing that happened removal must be programmed. In any case, we recommend the removal of the device before the recovery in the following cases: painful reaction due to the material, infection, implant rupture. Failure to carry out the removal and execution of a corrective action may result in complications.

### EXPOSURE TO THE ENVIRONMENT

Store in a dry place away from heat sources.

### DISPOSAL

The device must be disposed of as waste Hospital ..

MANUFACTURER

REFERENCE

BATCH NUMBER

WARNING:  
READ INSTRUCTION FOR USE

EXPIRATION DATE

SINGLE USE

DO NOT RESTERILIZE

KEEP AWAY FROM  
SUNLIGHT

KEEP DRY

DO NOT USE IF PACKAGING  
IS DAMAGED