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IMPORTANT INFORMATION FOR THE USE OF THE ENDOLOG SYSTEM

DESCRIPTION OF THE DEVICE

A device for the correction of hallux valgus in the technical dislocation of the metatarsal head.

Description: titanium nail, size 44 - Ø 4,5 x 40 mm length, size 45 - Ø 4,5 x 45 mm length, size 46 - Ø 46,9 x 4,5 mm length. - Elastic nail Endolog Adapta - Package of 1 pcs. STERILE.

The device Endolog adds screws of diameter 3.66 mm, in lengths of 15, 17, 20 or 25 mm pack of 1 pcs. STERILE. Material: Titanium alloy Ti 6Al-4V ELI -ISO 5832-3 - ASTM F136.

TERMS OF USE

The implant must be used by surgeons practitioners who have received the necessary information about the implantation of this device.

INDICATIONS

The system Endolog should be temporary used : fixation of the metatarsal head respect to the distal bone segment during surgical interventions of the distal osteotomy in hallux valgus correction surgery in patients of adult age. The choice of the size of the components of the system is the responsibility of the surgeon in function of the anatomy of the patient.

PERFORMANCE

The mechanical characteristics of this system must allow the synthesis of the distal osteotomy until the total bone consolidation. The performance can be hindered by a poor use of the material or from a deficient control postoperative patient.

SIDE EFFECTS

Pain and abnormal sensations due to the material, infection, neurological complication involving an eventual paralysis, nonunion, death.

PRECAUTIONS FOR USE

All the components of the Endolog system are supplied STERILE. Ethylene oxide (EO) sterilization complies with UNI EN ISO 11135-1. Use only if the package is intact and not damaged. Do not re-sterilize. SINGLE USE.

Pre-operative

A good knowledge of the products and of the surgical technique are essential. The choice of using the device is assigned to a surgeon competent according to the surgical indications and contraindications and the patient's anatomy.

Intra-operative

The instructions for use of the device Endolog must be followed carefully on the surgical technique. It 'not recommended to associate our implant to other similar devices.

Post-operative

A number of precautions must be given to the patient by the surgeon: limitation of his physical activity, compliance with the instructions of the surgeon.

CONTROL

Visual control of the assembly and through x-ray ; Control of the joint mobility; Controlling the closing of the different pieces of locking.

CONTRAINDICATIONS

active infections, fevers or leukocytes, obesity, inability to stabilize the system, mental illness, allergy to a metal, disagreement of the patient. The evaluation of which is deferred to the expertise of the surgeon, who decides if they are to be considered relative or absolute.

NUCLEAR MAGNETIC RESONANCE (NMR)

The device has not been tested for compatibility with diagnostic equipment Nuclear Magnetic Resonance ((MRI).

PACKAGING

Each component of the ENDOLOG device is provided STERILE and DISPOSABLE: will be used directly in the operating room. Reuse may cause a re-infection or cross-infection, or impair the functional performance of the device. Sterility is guaranteed only if the packaging is intact and not damaged. Do not re-sterilize.

ASSEMBLY OF THE DEVICE

The assembly of the Endolog system is made by the surgeon following the surgical technique, recommended in brochure illustration.

RADIOACTIVE DEVICES

The system Endolog is not radioactive.

CHANGE DEVICE PERFORMANCE

This system is expected to synthesize a distal metatarsal osteotomy until the bone healing. After this period, must be programmed ablation. In any case it is recommended that the ablation of the device before bone consolidation in the following cases: painful reaction due to the material, infection, rupture of a implant.

EXPOSURE TO THE ENVIRONMENT

Store in a dry place, away from sources of heat.

DISPOSAL

The device must be disposed of as Hospital waste.

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