LARK'S

SURGICAL DEVICES

dispositivo medico fabbricato in Italia
medical device manufactured in Italy

LARK'S

Endolog Adapta
HALLUX VALGUS ADAPTIVE IMPLANT

SURGICAL DEVICES SURGICAL DEVICES

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With over 10 years consolidated experience, the Endolog implant improves into a new system: the [Endolog Adapta](#).

The new product allows to further simplify the mini-invasive surgical treatment of hallux valgus deformity, making the implant adaptable to the metatarsal canal, ensuring shorter surgical time and, at the same time, excellent stability.

The [Endolog Adapta](#) implant composed of an elastic component with a curvilinear shape (intramedullary nail), which is produced in a single size and is adaptable to the metatarsal canal diameter, and an angle-stable screw.

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The implant is composed of an elastic curvilinear component manufactured in one size and an angle-stable screw.

The main feature of the system is the elastic curvilinear component which when introduced into the metatarsal canal after the metatarsal osteotomy, fits to the diameter of the metatarsal canal and achieves apposition to three distinct points:

- on the lateral cortex, with regard to the apex of the curvature,
- near the proximal metaphysis, with regard to the most proximal part of the implant,
- the lateral surface of metatarsal head, at the distal plate-shaped part of the implant.

The angular stability screw, introduced in the metatarsal head through the hole of the plate of the nail is the final element locking in stability of the implant.

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Endolog Adapta implant is proposed for the mini-invasive surgical treatment of hallux valgus deformity and for the correction of the follow entities:

- Media / Mild: IMA < 11°; MTF < 20°
- Moderata / Moderate: IMA > 11° < 16°; MTF > 20 < 40°
- Severa / Severe: IMA > 16° < 18°; MTF > 40°
- Revisione di precedenti interventi / Revision of previous surgery

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Due to the characteristics previously listed, the main feature of Endolog Adapta is the elastic stabilisation achieved between the resected metatarsal head and the 1st metatarsal. This leads to periosteal callus formation, allowing early mobilisation of the MTP joint.

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The implant provides a stable but not rigid distal osteotomy with significant improvements in all clinical parameters assessed in the pre-operative.

The implant also ensures a rapid normalisation post-operatively with:

- no pain
- good alignment
- disappearance of the areas of plantar hyperkeratosis
- no joint stiffness

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
The pre-operative planning should be assessed on radiographs in the upright position and assessing the following parameters:

- *Inter-metatarsal angle (IMA)*
- *Metatarsal Phalangeal angle (MTF)*
- *Proximal articular surface angle (PASA)*

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One very important step concerns the exostosectomy, which must be performed minimally as this is going to affect the degree of correction of the inter-metatarsal angle.

Exostosectomy must be made by taking distal-proximal wedge shape with a sagittal saw, in order to allow the implant to correct the valgus of the distal phalanx, thereby avoiding an additional surgical procedure.

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This wedge-shaped exostosectomy allows the correction of the PASA. (fig. 1)


Fig. 1

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If during the exostosectomy it is necessary to rotate the metatarsal head compared to the front section clockwise or counter clockwise, then the sagittal saw blade is angled appropriately oblique on the dorsal-plantar direction (fig.2).

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*The patient is under local regional anesthesia with tourniquet at level of foot neck.
In the case of severe valgus it is recommended to make adductor tenotomy performed with a micro scalpel blade.*

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A dorsal-medial, longitudinal incision of 2,0 / 2,5 cm is made corresponding to the exostosis of the 1st Metatarsal.



*The neuro-vascular bundle is isolated and protected.
A longitudinal incision is made to the MTP joint capsule.*

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
Disto-proximal exostosectomy made by taking a wedge shape from 2-4 mm to 0 mm.

The purpose of this important step allows you to restore the valgus of the distal phalanx in addition PASA once the flat, medial surface of the metatarsal head is resting on the flat part of the implant (see "pre-operative planning").

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A linear transverse osteotomy is performed at the level of the distal column of the first metatarsal according to pre-operative planning. An oblique osteotomy can be used to lengthen or shorten the metatarsal.

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Preparation of the access to the metatarsal canal with the size 46 test nail mounted on the impactor. This instrument, in addition to preparing the access of the implant into the canal, makes it possible to verify the efficacy of the displacement between the metatarsal head and 1st metatarsal of the implant will perform once it has been inserted into the metatarsal canal.

During this step, if it is necessary to reduce the displacement of the metatarsal head in relation to the first metatarsal bone, very gradually increase the exostectomy, with very small cut thicknesses, to achieve optimal displacement.

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*Assembling **Endolog Adapta** implant on the dedicated impactor, locking it with the appropriate screw.*

Care should be taken when inserting the implant to ensure it is without rotation and perpendicular to the longitudinal axis of the metatarsal.

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During insertion of the implant, displacement of the metatarsal head will occur with respect to the metatarsal. Continue with the insertion of the implant until the locking screw head is aligned with the centre of the metatarsal head.

Check the metatarsal head is in the correct position with fluoroscopy.

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
Note

With the implant in the correct position, manually adjust:

- the plantar position of the metatarsal head;
- **ensure the flat, medial surface of the metatarsal head is supported by the flat part of the implant. Never adjust the implant to fit against the metatarsal head.**

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WARNING:


The placement of the implant must be as vertical as possible to the longitudinal axis of the metatarsus, preventing rotation of the same during the placement.

In order to facilitate the early placement of the implant into the metatarsal canal will need to expand the osteotomy cut longitudinally exercising manual traction of the distal phalanx.

To complete implant placement with a hammer and by light blows sequences.

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Once these two manual adjustments have been made the metatarsal head can be temporarily stabilised with one diameter 1.0mm k-wire, inserted through the micro hole located on the implant, under the locking screw.

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
Remove the impactor from the implant by unscrewing the locking screw insert the angle-stable screw. The screw is self-drilling and self-tapping, it is not required to make the pilot hole.

During this step it is important to keep the parallax of the screw respect the implant, to allow proper grip to the implant.

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
Before closing the capsule, smooth the resected angle of the metatarsal to prevent soft-tissue irritation.

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Post-operative shoe with heel support for off-loading of the forefoot.

Begin with delicate movements of flexion/extension of the MTP joint, gradually increasing the rehabilitation.

First x-ray at 60 days.

After bone consolidation is complete the implant can be removed.

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Titanium Ti6Al4V ELI - ASTM F 136 - ISO 5838 - 5832-3

Surface treatment by anodic oxidation Type II

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Sterility is guaranteed only if the package is undamaged.
Do not re-sterilize. SINGLE USE ONLY.

BIBLIOGRAFIA / BIBLIOGRAPHY

J Orthopaed Traumatol (2009) 10 (Suppl 1):S41-S60

Hallux valgus mini-invasive surgery: comparison of two techniques - R.Valentini, G. Fancellu
Orthopaedic and Traumatologic Clinic, Trieste University (Trieste, IT)

8th EFAS Congress - 2nd-4th September 2010 - CMU, Geneva, Switzerland

A comparison between Endolog technique and Percutaneous Distal Metatarsal Osteotomy (PDO)
Carlo Biz, Roberto Aldegheri
Padua University - Orthopedic Clinic (Padua, IT)

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J Orthop Surg Res - 2015 Jul 2;10:102. doi: 10.1186/s13018-015-0245-1.

Endolog technique for correction of hallux valgus: a prospective study of 30 patients with 4-year follow-up
Carlo Biz, Marco Corradin, Ilaria Petretta, Roberto Aldegheri
Padua University - Orthopedic Clinic (Padua, IT)

Int Orthop - 2019 Mar;43(3):625-637. doi: 10.1007/s00264-018-4138-x. Epub 2018 Sep 14.

Minimally invasive surgery for hallux valgus: a systematic review of current surgical techniques
Francesc Malagelada, Cyrus Sahirad, Miki Dalmau-Pastor, Jordi Vega, Rej Bhumbra, María Cristina Manzanares-Céspedes
Royal London Hospital, UK

Foot Ankle Surg - 2021 Jan;27(1):46-51. doi: 10.1016/j.fas.2020.02.001. Epub 2020 Feb 5.

The Endolog technique for moderate to severe hallux valgus treatment: Clinical and radiographic analysis of 194 patients
Federico Bertolo, Andrea Pautasso, Claudio Cuocolo, Daniela Invernizzi, Francesco Atzori
Turin University - Orthopedic Clinic (Turin, IT)

Foot Ankle Int . 2021 Apr;42(4):409-424. doi: 10.1177/1071100720969676. Epub 2020 Dec 15.

Functional and Radiographic Outcomes of Minimally Invasive Intramedullary Nail Device (MIIND) for Moderate to Severe Hallux Valgus
Carlo Biz, Alberto Crimi, Ilaria Fantoni, Jacopo Tagliapietra, Pietro Ruggieri
Padua University - Orthopedic Clinic (Padua, IT)

Produttore e distributore / Manufacturer and distributor:

Medical Due S.r.l.
Via Trento, 43
IT-37014 Castelnuovo del Garda (VR) - Italy
Tel. +39 (0) 45 6450777 - fax +39 (0) 45 8525203
www.medicaldue.it
P.IVA IT-02895150239 - Reg.Imp. VR1998-63861
R.E.A. 292702 - Cap. soc. Euro 10.400,00 i.v.