**LARK'S**  
SURGICAL DEVICES

LARK'S

**KONE-LOG System**

DISPOSITIVO PER ARTRORISI ENDOSENOTARSICA  
*SINUS TARSII ARTHROEREISIS DEVICE*

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*The KONE-LOG device was developed to correct and restrict excessive pronation of flatfoot in adults and children.*

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*Conical shape with not traumatic thread.*

*Diameters of between 6.0 and 13 mm*

*Titanium alloy with Type II anodizing surface treatment.*

*Cannulation 2.0 mm*


*The devices are provided in sterile individual packs and are intended for single use.*

*The instrument set includes a guide wire, sinus tarsi spacer (a), the set of cannulated test implant dilators with diameters of 7.0 to 13 mm (b), a cannulated screwdriver for implant insertion (c) and a removal screwdriver (d).*

*Attention: the surgeon must have adequate training and an in-depth familiarity with use of the Kone-Log implant and instrument set.*

*Knowledge of the surgical technique is extremely important.*

*Preoperative planning is fundamental.*

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*The use of Kone-Log device is indicated for the treatment of hyperpronation of the foot and stabilisation of the subtalar joint. The system is designed to prevent forward, downward or medial displacement of the talus, in order to recover normal subtalar joint mobility, whilst preventing excessive pronation and the corresponding consequences.*

*The indications for use include:*

- *severe pronation;*
- *calcaneal valgus deformity;*
- *plantar flexed talus;*
- *failed correction with long-term orthotic treatment;*
- *painful and congenital flatfoot;*
- *post-tarsal coalition repair;*
- *subtalar instability;*
- *posterior tibial tendon dysfunction;*
- *paralytic flatfoot deformity.*

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*Use of the Kone-Log device is contraindicated in the following cases:*

- *Active or suspected sepsis or in immunocompromised patients*
- *In patients with prior sensitisation to titanium;*
- *In patients with inadequate bone stock;*
- *In patients with superstructural alignment deformity or in patients with certain metabolic disorders;*
- *Patient disagreement.*

*The evaluation of the above aspects is the responsibility of the surgeon, who must decide whether they are relative or absolute contraindications.*

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*Perform a linear incision approximately 2-3 cm long above the tarsal sinus on the lateral aspect of the foot. Identify the intermediate dorsal cutaneous nerve and retract with care.*

- *Perform a blunt dissection along the fascia and the capsule covering the tarsal sinus and make an incision.*

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*For easy inserting of the guide wire, to increase the space between the heel and astragalus with a special instrument.*

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*Insert the guide wire into the tarsal canal, until the end of it comes into contact with the soft tissue on the medial part of the foot.*

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*Introduce the cannulated test implant dilators in sequence, starting from that with the smallest diameter, until the desired correction is obtained.*

- *Check the intraoperative position by fluoroscopy.*
- *Having selected the correct dilator diameter, which will be identical to the implant to be used, remove the dilator, leaving the guide wire in position.*
- *Choose the implant diameter that matches that of the last dilator used.*

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
*Place the device on the guide wire and, using the cannulated hexagonal screwdriver, insert to the desired depth inside the subtalar joint.*


- *Check the necessary correction has been achieved by moving the subtalar joint and confirm it by intraoperative x-ray or fluoroscopic imaging in both anteroposterior and laterolateral projections.*
- *Remove the cannulated screwdriver and guide wire.*
- *Analyse the movement of the subtalar joint.*
- *Irrigate thoroughly with sterile saline solution.*
- *Suture the capsule and deep and superficial tissues according to standard procedure.*

**FIG.4**

**FIG.5**

## POST OPERATORIO / POST OPERATIVE-CARE

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*Assuming no adjunctive procedures were performed, a protective, weight bearing, below the knee walking cast or boot for 2-4 weeks is used. A gradual return to limited activity in 4-6 weeks is permitted as tolerated.*

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

When the screwdriver meets resistance, continue turning the screwdriver (**FIG. 6**), pulling it slowly until the device has been completely removed (**FIG. 7**).

- Irrigate thoroughly with sterile saline solution.
- Suture the capsule and deep and superficial tissues according to standard procedure.

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*Sterility is guaranteed only if the package is undamaged.*

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**Fabbricante e distributore**  
***Manufacturer and distributor***

Medical Due S.r.l.  
Via Trento, 43  
IT-37014 Castelnuovo del Garda (VR) - Italy  
Tel. +39 (0) 45 6450777 - fax +39 (0) 45 8525203  
[www.medicaldue.it](http://www.medicaldue.it)  
P.IVA IT-02895150239 - Reg.Imp. VR1998-63861  
R.E.A. 292702 - Cap. soc. Euro 10.400,00 i.v.

Rev 01\_2021