

LARK'S

SURGICAL DEVICES

LARK'S

HAND-LOG

Immobilizzatore per la chirurgia della mano
Supporto malleabile rivestito in silicone medical grade



Cod: H-L

- Senza piombo e Latex-free.
- Radiotrasparente - idoneo per l'utilizzo con il fluoroscopio.
- Previene la rotazione della mano e del polso durante l'intervento.
- Con ala anti ribaltamento.
- Dopo l'intervento, facile da riposizionare nella sua forma originale.
- Lamina di alluminio malleabile rivestita in silicone medical grade.
- Facile posizionamento della mano e comfort per il paziente.
- Riutilizzabile.
- Autoclavabile.
- Universale, ambidestro.
- Conforme alle direttive EC.



Produttore e distributore

Medical Due S.r.l. - Via Trento, 43 - 37014 Castelnuovo del Garda (VR)
Tel. 045 6450777 - fax 045 6450798 - www.medicaldue.it

P.IVA IT-02895150239 - Reg.Imp. VR1998-63861 - R.E.A. 292702 - Cap. soc. Euro 10.400,00 i.v.

Rev 01_2011

LARK'S

SURGICAL DEVICES

LARK'S

HAND-LOG

Immobilizer for hand surgery
Pliabe sheet coated of silicone medical grade



CE

Cod: H-L

- Lead-free and Latex-free.
- Radiolucent - Allows for the use of fluoroscopy.
- Prevents arm and hand rotation during surgery.
- With wing anti-rollover.
- After surgery, easily repositionable in its original form.
- Pliable aluminium sheet coated of silicone medical grade.
- Easy hand positioning and for patient comfort.
- Reusable.
- Autoclavable.
- Universal, right and left.
- Complies with EC Directives.



Manufacturer and distributor

Medical Due S.r.l. - Via Trento, 43 - IT 37014 Castelnuovo del Garda (VR) - Italy
Tel. +39 (0) 45 6450777 - fax +39 (0) 45 6450798 - www.medicaldue.it
P.IVA IT-02895150239 - Reg.Imp. VR1998-63861 - R.E.A. 292702 - Cap. soc. Euro 10.400,00 i.v.

Rev 01_2011